

REMP Installations

First -80°C BioBank

This new automated storage and retrieval system from REMP is a key component of Pfizer's BioBank and is a result of customer focused collaboration leading to state-of-the-art storage solutions.

Utilizing the enormous potential of biological samples

Thanks to the sequencing of the human genome and continually improving analytical methods in the field of molecular biology, biological samples have acquired a new significance and value. Where samples were once gathered for specifically stated purposes, had limited use, and were stored for archival purposes, they can now be stored under the best possible conditions, for multiple uses in further studies.

Modern investigative methods in the fields of "personalized medicine", pharmacogenomics, pharmacometabonomics and proteomics now provide far better insight into the modes of action of different therapeutics, enabling them to be optimized and administered more easily and accurately. Properly stored clinical trial samples have the potential to expedite development of investigational

new medicines, thereby accelerating their delivery to patients. The biological samples can help researchers uncover pathways and metabolic processes which may have been previously inaccessible or undiscovered.

The pharmaceutical industry has recognized this enormous research potential and the value to both scientists and patients on making these study samples available worldwide. This requires samples to be stored under precisely defined conditions at extremely low temperatures. They should not be subjected to freeze/thaw cycles and any temperature fluctuations within the system must be minimized. The multitude of samples and defined storage conditions practically rule out any kind of manual handling of samples. This critical storage process needs to be automated in order to obtain samples of the best possible quality.

Pfizer's Big Chill

Pfizer is the world's largest pharmaceutical company with the largest privately-financed research organization. Over 115,000 employees in more than 150 countries work towards improving health. Its business activities encompass the fields of prescription medicines, over-the-counter medicines and veterinary medicines. Thanks to continuous research and development as well as long-term planning, Pfizer is repeatedly successful in being able to offer effective medicines in the fight against and alleviation of diseases.

Pfizer also makes an important contribution to the provision of Switzerland's healthcare with over 150 medicines.



Fully automated storage at -80°C and -20°C

REMP has had long-standing experience in the design and construction of storage systems for chemical compounds.

The objective of the Pfizer BioBank project was to ensure the precisely defined, fully automated storage of biological samples at temperatures of -80°C and -20°C , and to safeguard the best possible quality of such samples over extended periods of time.



Photos: Tim Cook/The Day Publishing Company (New London, CT)



Based on proven storage technology

The first REMP Bio-Sample Store™ was developed jointly with Pfizer Global Research and Development, the R&D arm of Pfizer, Inc., and is able to store millions of samples at -80°C and -20°C, concurrently, and under controlled dry air conditions to prevent frost build up. For the design, the already proven concept of REMP's -20°C storage solutions for chemical compounds was modified by fitting additional large, automated freezer compartments that operate at -80°C. New door mechanisms to ensure the time-critical opening of the compartments were also developed. In order to avoid mechanical problems, due to the low temperature environment, only passive mechanical components were used. Newly developed software algorithms further guarantee optimally selected storage locations and best possible

access to samples. The integrated REMP Tube Punching Module™ ensures that cherry-picking or reformatting of samples always occurs at -20°C, within the Bio-Sample Store. The storage unit has fully redundant cooling and robotic systems, which allows service and maintenance to occur without any loss of system use or effect on the controlled storage environment.

The world's first large scale automated -80°C storage facility

Pfizer's BioBank is the world's first large scale automated store for biological samples in which samples can be stored at both -80°C and -20°C. Equipped with REMP's patented Tube Technology™, Pfizer operates a complete system that minimizes the number of freeze/thaw cycles subjected to a sample, thereby maximizing optimal sample quality. All of the basic conditions were developed

in cooperation with Pfizer according to their requirements. Thanks to its long-standing experience in building automated storage facilities and managing large scale projects, REMP was the first company to setup an automated -80°C store in record time. The store was installed and commissioned on schedule.

Key dates:

- Contract signed – June 2004;
- FAT (factory acceptance test) – April 2005;
- SAT (site acceptance test) – March 2006.

REMP AG

Weststrasse 12
CH-3672 Oberdiessbach
Switzerland
Tel.: +41 31 770 70 70
www.remp.com

Tecan US Inc.

REMP Business Unit
150 Hopping Brook Road
Holliston, MA 01746
USA
Tel.: +1 508 429 2200

Tecan Japan Co., Ltd.

REMP Business Unit
Kawasaki-Tech Center Bldg. 17F
580-16, Horikawa-cho
Saiwai-ku, Kawasaki-shi 212-0013
Tel.: +81 44 542 7021

REMP® is a registered trademark of REMP AG, Oberdiessbach Switzerland. REMP Bio-Sample Store™, REMP Tube Punching Module™ and REMP Tube Technology™ are trademarks of REMP AG, Oberdiessbach, Switzerland.